

BeneCheck™

Uric Acid Monitoring System



BeneCheck's Electrochemical Uric Acid Monitoring System is the best choice for self-diagnostic.

Accurate

Easier to use

Faster results

BeneCheck Uric Acid monitoring system is brought to you by General Life Biotechnology Co., Ltd. The BeneCheck System uses leading amperometric technology in a unique electrochemical process. Capillary sampling method and sample zone window design, provides you user friendly operation. The fast, accurate and powerful system helps you to monitor your uric acid level easily.


BeneCheck™

Uric Acid Monitoring System

SPECIFICATION

| | | | |
|-------------------------|-------------------------------------|-------------------|------------------------------------|
| Principle | Electrochemical Biosensor | Relative Humidity | Less than 95% |
| Test Sample | Capillary Whole Blood | Memory | 180 test results |
| Measuring Time | ≤ 15 seconds | Battery Life | Approximately 1,000 tests |
| Measuring Range | 3 - 20 mg/dL (0.18 - 1.19 mmol/L) | Battery Type | 3V (CR2032) Lithium Battery x 1 |
| Sample Volume | ≥ 0.7 μL | Dimensions | 61*77*19 mm |
| Strip Storage Condition | 10 - 30°C | Weight | About 47g (battery not included) |
| Operation Temperature | 10 - 40°C | | |

PACKAGING CONTENTS




Uric Acid Meter



UA
Test Strips



Lancing Device



Carrying Bag




Lancet

STRIPS WARRANTY

Unopened Shelf life : 18 months
(From manufacturing Day)

Opened Shelf life : 3 months after first opening

We provide a momentum to improve human health and life quality for all



BeneCheck™

General Life Biotechnology Co., Ltd.
5F, No. 240, Shinshu Road,
Shin Juang City, Taipei, 242 Taiwan

TEL : +886-2-2205-3698

FAX : +886-2-2208-0779

E-mail: service@glbiotech.com.tw

www.glbiotech.com.tw
www.BeneCheck.com.tw

