

BeneCheck™

Blood Glucose Monitoring System



BeneCheck's Electrochemical Blood Glucose Monitoring System is the best choice for self-diagnostic.

Easier to use

Faster results

Painless procedure

BeneCheck Blood Glucose Monitoring System is brought to you by General Life Biotechnology Co., Ltd. The BeneCheck System use leading amperometric technology in an unque electrochemical process. Capillary sampling method and sample zone window design, provides you user friendly operation. Detection in a few pain. This fast, accurate and powerful system helps you to monitor your blood glucose level easily.

BeneCheck™

Blood Glucose Monitoring System


SPECIFICATION

Principle	Electrochemical Biosensor	Relative Humidity	Less than 95%
Test Sample	Capillary Whole Blood	Memory	180 test results
Measuring Time	≤ 10 seconds	Battery Life	Approximately 1,000 tests
Measuring Range	18 - 600 mg/dL (1.0 - 33.3 mmol/L)	Battery Type	3V (CR2032) Lithium Battery x 1
Sample Volume	≥ 0.7 μL	Dimensions	61*77*19 mm
Strip Storage Condition	10 - 30°C	Weight	About 47g (battery not included)
Operation Temperature	10 - 40°C		

PACKAGING CONTENTS



Glucose Meter



GLU
Test Strips



Lancing Device



Carrying Bag




Lancet

STRIPS WARRANTY

Unopened Shelf life : 24 months
(From manufacturing Day)

Opened Shelf life : 3 months after first opening

We provide a momentum to improve human health and life quality for all



BeneCheck™

General Life Biotechnology Co., Ltd.
5F, No. 240, Shinshu Road,
Shin Juang City, Taipei, 242 Taiwan

TEL : +886-2-2205-3698

FAX : +886-2-2208-0779

E-mail: service@glbiotech.com.tw

www.glbiotech.com.tw

www.BeneCheck.com.tw

